

200 µL, Racked, Low Retention, Sterile, 96 tips/rack

Specification Sheet (Part/REF # 68008-22)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

Product Specifications

Disposable Tip Volume	200 µL
Tip Length	Standard
Packaging	Racked
Quantity Per Each	960
Filtered	No
Low Retention	Yes
Sterile	Yes
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recomended Pipette Volume Ranges	20 µL, 25 µL, 50 µL, 100 µL, and 200 µL Pipettes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68008-22>

Spec. Version B