

300 μ L, Racked, Low Retention, 96 tips/rack

Specification Sheet (Part/REF # 68009-21)

Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, pre-sterilized, and macro volume.

Product Specifications

Disposable Tip Volume	300 μ L
Tip Length	Standard
Packaging	Racked
Quantity Per Each	960
Filtered	No
Low Retention	Yes
Pre-Sterilized	No
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/pipettes/68009-21>