

300 µL, Racked, Filter, Sterile, 96 tips/rack

Specification Sheet (Part/REF # 68009-13)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

Product Specifications

Disposable Tip Volume	300 µL
Tip Length	Standard
Packaging	Racked
Quantity Per Each	960
Filtered	Yes
Low Retention	No
Sterile	Yes
Recomended Pipette Volume Ranges	250 µL and 300 µL Pipettes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/pipettes/68009-13>

Spec. Version B