

# PB600-1 Dispenser

Specification Sheet (Part/REF # 83700)



The VDRL Antigen Syringe is designed to simplify the dispensing steps for VDRL antigen testing. The test is a microscopic nontreponemal flocculation test to be used for the detection of reagin.

## Product Specifications

**Compatible with** 700 and 1700 series syringes in the following volumes: 10  $\mu$ L, 25  $\mu$ L, 50  $\mu$ L, 100  $\mu$ L, 250  $\mu$ L, 500  $\mu$ L, 1 mL, and 2.5 mL

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/other-laboratory/83700>

**Spec. Version B**