

**2.5 mL, Model 1002 LTN CTC SYR,
Cemented NDL, 23 ga, point style 5
(110 °C)**

Specification Sheet (Part/REF # 203084)



Product Specifications

Volume	2.5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Needle (LTN)
Needle	Cemented Needle
Gauge	23 gauge
Point Style	5: Conical Side Port
Needle Length	56 mm (2.20 inches)
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	29 psig (2bar)
Maximum Temperature	230°F (110°C)
Minimum Temperature	68°F (20°C)
Needle Inner Diameter	0.008 inches (0.20 mm)
Needle Outer Diameter	0.025 inches (0.64 mm)
Std Needle Dead Volume	4.53 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/syringes/203084>

Spec. Version B