

HVDP 2-5 Valve Black with Female Luer

Specification Sheet (Part/REF # 86789)



These Hamilton Gastight syringes are ideal for Purge and Trap applications because you can easily remove the syringe plunger and load water samples into the syringe barrel. You can add an internal standard or surrogate sample directly to your sample through the twist valve with a 10 µL syringe.

Product Specifications

Sales Quantity	1 valve
----------------	---------

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/syringes/86789>

Spec. Version B