

10 mL, Model 1010 RN SYR, Needle Sold Separately

Specification Sheet (Part/REF # 7652-01)



The 1000 series syringes are a mid-volume Gastight syringe. This series of syringes is ideal for dispensing volumes from 1 mL up to 100 mL. A variety of different terminations and needle options are available in this series.

Product Specifications

| | |
|------------------------------|---|
| Volume | 10 mL |
| Syringe Type | Gastight |
| Syringe Series | 1000 |
| Termination | Removable Needle (RN) |
| Needle | Large RN |
| Gauge | No Needle Included |
| Point Style | No Needle Included |
| Needle Length | No Needle Included |
| Fluid Path | Borosilicate glass, PTFE SST, Borosilicate glass, PTFE |
| Autoclavable | Yes |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 115°C / 239°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.574 inches (14.6 mm) |

| | |
|-------------------------------|------------------------|
| Barrel Outer Diameter | 0.695 inches (17.7 mm) |
| Std Needle Dead Volume | N/A |
| Sales Quantity | 1 syringe per pack |

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/syringes/7652-01>