

250 µL, Model 1725 AD-SAL Syringe, M8 Threads

Specification Sheet (Part/REF # 203310)



The 1700 series syringes are the gas tight version of the original 700 series syringes. This series of gas tight syringes is ideal for dispensing volumes from 1 µL up to 500 µL. A variety of different terminations and needle options are available in this series.

Product Specifications

Needle Sold Separately	Yes
Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	AccuDil® (AD)
Needle	No Needle Available
Plunger Style	Saltline (SAL)
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Outer Diameter	0.028 inch (0.7 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/syringes/203310>

Spec. Version B