

Custom Length, Heavy Wall PTFE Tubing (0.123" O.D. / 0.063" I.D.)

Specification Sheet (Part/REF # 20942)



Product Specifications

Tubing Material	PTFE
Tubing Inner Diameter (I.D.)	0.063 in / 1.600 mm
Tubing Outer Diameter (O.D.)	0.123 in / 3.124 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/20942>

Spec. Version B