

Custom Length PTFE Tubing, 12 Gauge (0.110" O.D. / 0.086" I.D.)

Specification Sheet (Part/REF # 20912)



Product Specifications

Tubing Material	PTFE
Gauge	12 gauge
Tubing Inner Diameter (I.D.)	0.086 in / 2.184 mm
Tubing Outer Diameter (O.D.)	0.110 in / 2.794 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/20912>

Spec. Version B