

# Custom Length PTFE Tubing, Heavy Wall, 19 gauge (0.063" O.D. / 0.031" I.D.)

Specification Sheet (Part/REF # 20943)



## Product Specifications

Tubing Material	PTFE
Gauge	19 gauge
Tubing Inner Diameter (I.D.)	0.031 in / 0.787 mm
Tubing Outer Diameter (O.D.)	0.063 in / 1.600 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/20943>

Spec. Version B