

Custom Length PTFE Tubing, 17 Gauge (0.073" O.D. / 0.049" I.D.)

Specification Sheet (Part/REF # 20917)



Product Specifications

Tubing Material	PTFE
Gauge	17 gauge
Tubing Inner Diameter (I.D.)	0.049 in / 1.245 mm
Tubing Outer Diameter (O.D.)	0.073 in / 1.854 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/20917>

Spec. Version B