

# Custom Length PTFE Tubing (0.155" O.D. / 0.125" I.D.)

Specification Sheet (Part/REF # 20941)



## Product Specifications

Tubing Material	PTFE
Tubing Inner Diameter (I.D.)	0.125 in / 3.175 mm
Tubing Outer Diameter (O.D.)	0.155 in / 3.937 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/20941>

Spec. Version B