

# Custom Length PTFE Tubing, 20 Gauge (0.059" O.D. / 0.035" I.D.)

Specification Sheet (Part/REF # 20920)



## Product Specifications

Tubing Material	PTFE
Gauge	20 gauge
Tubing Inner Diameter (I.D.)	0.035 in / 0.889 mm
Tubing Outer Diameter (O.D.)	0.059 in / 1.499 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/20920>

Spec. Version B