

# Custom Length PTFE Tubing, 24 Gauge (0.043" O.D. / 0.023" I.D.)

Specification Sheet (Part/REF # 20924)



## Product Specifications

Tubing Material	PTFE
Gauge	24 gauge
Tubing Inner Diameter (I.D.)	0.023 in / 0.584 mm
Tubing Outer Diameter (O.D.)	0.043 in / 1.092 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/20924>

Spec. Version B