

Custom Length PTFE Tubing, 26 Gauge (0.037" O.D. / 0.019" I.D.)

Specification Sheet (Part/REF # 20926)



Product Specifications

| | |
|------------------------------|---------------------|
| Tubing Material | PTFE |
| Gauge | 26 gauge |
| Tubing Inner Diameter (I.D.) | 0.019 in / 0.483 mm |
| Tubing Outer Diameter (O.D.) | 0.037 in / 0.940 mm |
| Tubing Length | Custom Length |
| Fitting Material | No Hub |
| Sales Quantity | 1 each |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/20926>

Spec. Version B