

# Custom Length PTFE Tubing, 30 Gauge (0.031" O.D. / 0.013" I.D.)

Specification Sheet (Part/REF # 20930)



## Product Specifications

Tubing Material	PTFE
Gauge	30 gauge
Tubing Inner Diameter (I.D.)	0.013 in / 0.330 mm
Tubing Outer Diameter (O.D.)	0.031 in / 0.787 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/20930>

Spec. Version B