

# M6 Fitting (0.134" I.D.) 1/pk

Specification Sheet (Part/REF # 32857)

For use with 1/8" O.D. Flanged Tubing

## Product Specifications

Fitting Material	CTFE
Fitting Inner Diameter (I.D.)	0.134 in / 3.404 mm
Fitting Outer Diameter (O.D.)	0.173 in / 4.394 mm
Threads	M6
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/32857>

Spec. Version B