

# Custom Length FEP Tubing, 18 Gauge (0.075" O.D. / 0.043" I.D.)

Specification Sheet (Part/REF # 20937)



## Product Specifications

Tubing Material	FEP
Gauge	18 gauge
Tubing Inner Diameter (I.D.)	0.043 in / 1.092 mm
Tubing Outer Diameter (O.D.)	0.075 in / 1.905 mm
Tubing Length	Custom Length
Fitting Material	No Hub

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/ fittings/20937>

Spec. Version B