

# Custom Length FEP Tubing, 22 Gauge (0.074" O.D. / 0.028" I.D.)

Specification Sheet (Part/REF # 20939)



## Product Specifications

Tubing Material	FEP
Gauge	22 gauge
Tubing Inner Diameter (I.D.)	0.028 in / 0.711 mm
Tubing Outer Diameter (O.D.)	0.074 in / 1.880 mm
Tubing Length	Custom Length
Fitting Material	No Hub
Maximum Temperature	200°C / 392°F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/ fittings/20939>

Spec. Version B