

Female Luer Lock, One Fitting with 12" PTFE Tubing, 13 Gauge (0.101" O.D. / 0.077" I.D.)

Specification Sheet (Part/REF # 90613)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	13 gauge
Tubing Inner Diameter (I.D.)	0.077 in / 1.956 mm
Tubing Outer Diameter (O.D.)	0.101 in / 2.565 mm
Tubing Length	12 inches
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90613>

Spec. Version B