

Female Luer Lock, One Fitting with 12" PTFE Tubing, 24 Gauge (0.043" O.D. / 0.023" I.D.)

Specification Sheet (Part/REF # 90624)



Product Specifications

| | |
|------------------------------|---------------------|
| Fitting Material | CTFE |
| Tubing Material | PTFE |
| Gauge | 24 gauge |
| Tubing Inner Diameter (I.D.) | 0.023 in / 0.584 mm |
| Tubing Outer Diameter (O.D.) | 0.043 in / 1.092 mm |
| Tubing Length | 12 inches |
| Number of Fittings | 1 |
| Fitting Type | Female Luer |
| Sales Quantity | 1 each |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90624>

Spec. Version B