

**Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 11 Gauge (0.120" O.D. / 0.096" I.D.)**

Specification Sheet (Part/REF # 90661)



**Product Specifications**

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	11 gauge
Tubing Inner Diameter (I.D.)	0.096 in / 2.438 mm
Tubing Outer Diameter (O.D.)	0.120 in / 3.048 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90661>

**Spec. Version B**