

# Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 13 Gauge (0.101" O.D. / 0.077" I.D.)

Specification Sheet (Part/REF # 90663)



## Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	13 gauge
Tubing Inner Diameter (I.D.)	0.077 in / 1.956 mm
Tubing Outer Diameter (O.D.)	0.101 in / 2.565 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90663>

Spec. Version B