

Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 14 Gauge (0.092" O.D. / 0.068" I.D.)

Specification Sheet (Part/REF # 90664)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	14 gauge
Tubing Inner Diameter (I.D.)	0.068 in / 1.727 mm
Tubing Outer Diameter (O.D.)	0.092 in / 2.337 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90664>

Spec. Version B