

# Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 15 Gauge (0.085" O.D. / 0.061" I.D.)

Specification Sheet (Part/REF # 90665)



## Product Specifications

|                              |                     |
|------------------------------|---------------------|
| Fitting Material             | CTFE                |
| Tubing Material              | PTFE                |
| Gauge                        | 15 gauge            |
| Tubing Inner Diameter (I.D.) | 0.061 in / 1.549 mm |
| Tubing Outer Diameter (O.D.) | 0.085 in / 2.159 mm |
| Tubing Length                | Custom Length       |
| Number of Fittings           | 1                   |
| Fitting Type                 | Female Luer         |
| Sales Quantity               | 1 each              |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90665>

Spec. Version B