

**Female Luer Lock, One Fitting with  
Custom Length PTFE Tubing, 16  
Gauge (0.079" O.D. / 0.055" I.D.)**

Specification Sheet (Part/REF # 90666)



**Product Specifications**

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	16 gauge
Tubing Inner Diameter (I.D.)	0.055 in / 1.397 mm
Tubing Outer Diameter (O.D.)	0.079 in / 2.007 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90666>

**Spec. Version B**