

# Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 18 Gauge (0.067" O.D. / 0.043" I.D.)

Specification Sheet (Part/REF # 90668)



## Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	18 gauge
Tubing Inner Diameter (I.D.)	0.043 in / 1.092 mm
Tubing Outer Diameter (O.D.)	0.067 in / 1.702 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90668>

Spec. Version B