

# One Hub Assembled with Custom Length Tubing (1/8")

Specification Sheet (Part/REF # 90691)

## Product Specifications

<b>Fitting Material</b>	CTFE
<b>Tubing Material</b>	PTFE
<b>Gauge</b>	1/8"
<b>Length</b>	Custom Length
<b>Sales Quantity</b>	1 each
<b>Fitting Type</b>	Female Luer
<b>Number of Fittings</b>	one

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/fittings/90691>

**Spec. Version B**