

Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 19 Gauge (0.063" O.D. / 0.039" I.D.)

Specification Sheet (Part/REF # 90669)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	19 gauge
Tubing Inner Diameter (I.D.)	0.039 in / 0.991 mm
Tubing Outer Diameter (O.D.)	0.063 in / 1.600 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90669>

Spec. Version B