

Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 20 Gauge (0.059" O.D. / 0.035" I.D.)

Specification Sheet (Part/REF # 90670)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	20 gauge
Tubing Inner Diameter (I.D.)	0.035 in / 0.889 mm
Tubing Outer Diameter (O.D.)	0.059 in / 1.499 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90670>

Spec. Version B