

Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 22 Gauge (0.049" O.D. / 0.029" I.D.)

Specification Sheet (Part/REF # 90672)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	22 gauge
Tubing Inner Diameter (I.D.)	0.029 in / 0.737 mm
Tubing Outer Diameter (O.D.)	0.049 in / 1.245 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90672>

Spec. Version B