

Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 24 Gauge (0.043" O.D. / 0.023" I.D.)

Specification Sheet (Part/REF # 90674)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	24 gauge
Tubing Inner Diameter (I.D.)	0.023 in / 0.584 mm
Tubing Outer Diameter (O.D.)	0.043 in / 1.092 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90674>

Spec. Version B