

# Female Luer Lock, One Fitting with Custom Length PTFE Tubing, 26 Gauge (0.037" O.D. / 0.019" I.D.)

Specification Sheet (Part/REF # 90676)



## Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	26 gauge
Tubing Inner Diameter (I.D.)	0.019 in / 0.483 mm
Tubing Outer Diameter (O.D.)	0.037 in / 0.940 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90676>

Spec. Version B