

**Female Luer Lock Tubing Assembly,  
One Fitting with Custom Length  
PTFE Tubing, 30 Gauge (0.031" O.D.  
/ 0.013" I.D.)**

Specification Sheet (Part/REF # 90680)



Product Specifications

|                              |                     |
|------------------------------|---------------------|
| Fitting Material             | CTFE                |
| Tubing Material              | PTFE                |
| Gauge                        | 30 gauge            |
| Tubing Inner Diameter (I.D.) | 0.013 in / 0.330 mm |
| Tubing Outer Diameter (O.D.) | 0.031 in / 0.787 mm |
| Tubing Length                | Custom Length       |
| Number of Fittings           | 1                   |
| Fitting Type                 | Female Luer         |
| Sales Quantity               | 1 each              |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/fittings/90680>

Spec. Version B