

**Female Luer Lock Tubing Assembly,
One Fitting with Custom Length
PTFE Tubing, 7 Gauge (0.180" O.D. /
0.150" I.D.)**

Specification Sheet (Part/REF # 90657)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	7 gauge
Tubing Inner Diameter (I.D.)	0.150 in / 3.810 mm
Tubing Outer Diameter (O.D.)	0.180 in / 4.572 mm
Tubing Length	Custom Length
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/90657>

Spec. Version B