

Female Luer Lock Tubing Assembly, One Fitting with 12" PTFE Tubing, 8 Gauge (0.164" O.D. / 0.134" I.D.)

Specification Sheet (Part/REF # 90608)



Product Specifications

Fitting Material	CTFE
Tubing Material	PTFE
Gauge	8 gauge
Tubing Inner Diameter (I.D.)	0.134 in / 3.404 mm
Tubing Outer Diameter (O.D.)	0.164 in / 4.166 mm
Tubing Length	12 inches
Number of Fittings	1
Fitting Type	Female Luer
Sales Quantity	1 each

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/fittings/90608>

Spec. Version B