

# Two Hubs Assembled with Custom Length Tubing (1/16" O.D., 1/32" I.D.) IS/IS

Specification Sheet (Part/REF # 88902)

## Product Specifications

<b>Hub</b>	Kel-F
<b>Length</b>	Custom Length
<b>Sales Quantity</b>	1 each
<b>Tubing Material</b>	PTFE
<b>Number of Hubs</b>	2
<b>Thread Type</b>	1/4"-28
<b>Tubing Outer Diameter</b>	0.063"
<b>Tubing Inner Diameter</b>	0.031"
<b>Nut Length</b>	IS (0.60")

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/fittings/88902>