

Two Hubs Assembled with Custom Length PTFE Tubing (1/16"O.D./ 0.042 I.D.) IS/IS with Strain Relief

Specification Sheet (Part/REF # 88901)

Product Specifications

Hub	Kel-F
Length	Custom Length
Sales Quantity	1 each
Tubing Material	PTFE
Number of Hubs	2
Tubing Outer Diameter	0.063"
Tubing Inner Diameter	0.042"
Thread Type	1/4"-28
Nut Length	IS (0.60")

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/fittings/88901>