

1/4"-28 Tubing Assembly, Two Fittings with Custom Length PTFE Tubing, (0.067" O.D. / 0.043" I.D.), IS/IS

Specification Sheet (Part/REF # 88900)



Product Specifications

Tubing Material	PTFE
Washer Material	Stainless Steel
Gauge	18 gauge
Tubing Inner Diameter (I.D.)	0.043 in / 1.092 mm
Tubing Outer Diameter (O.D.)	0.067 in / 1.702 mm
Washer Inner Diameter (I.D.)	0.073 in / 1.854 mm
Washer Outer Diameter (O.D.)	0.189 in / 4.801 mm
Tubing Length	Custom Length
Nut Length	IS (0.60")
Fitting Material	CTFE
Number of Fittings	2
Fitting Inner Diameter (I.D.)	0.089 in / 2.261 mm
Fitting Outer Diameter (O.D.)	0.118 in / 2.997 mm
Threads	1/4"-28
Sales Quantity	1 each

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/fittings/88900>

Spec. Version B