

Two Hubs Assembled with Custom Length PTFE Tubing (O.D. 0.120" / I.D. 0.091") IS/IS

Specification Sheet (Part/REF # 88910)

Product Specifications

Hub	Kel-F
Length	Custom Length
Sales Quantity	1 each
Tubing Material	PTFE
Number of Hubs	2
Thread Type	1/4"-28
Tubing Outer Diameter	0.120"
Tubing Inner Diameter	0.091"
Nut Length	IS (0.60")

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/fittings/88910>