

# Two Hubs Assembled with Custom Length PTFE Tubing (1/8" O.D., 1/16" I.D.) IS/IS with Strain Relief

Specification Sheet (Part/REF # 88913)

## Product Specifications

Hub	Kel-F
Length	Custom Length
Sales Quantity	1 each
Tubing Material	PTFE
Number of Hubs	2
Thread Type	1/4"-28
Tubing Outer Diameter	0.125" O.D.
Tubing Inner Diameter	0.063" I.D.
Nut Length	IS (0.60")

Specifications are subject to change without notice

**Spec. Version B**

<https://www.hamiltoncompany.com/laboratory-products/fittings/88913>