

VALVE ASSY,DILUTER,ML560

Specification Sheet (Part/REF # 53175-01)

Product Specifications

Valve Model	ML500 Valves
Valve Model Description	Diluter, dual valve with slots for the ML560
Valve Size	DUAL
Valve Outer Body Dimension	DUAL
Actuated	ML560
Break Away Torque	15 oz-in (105.9 mNm)
Panel Mounting	No
Operating Temp High	38 °C (100 °F)
Operating Temp Low	-12 °C (10 °F)
Operating Pressure	0 to 100 psig (0 to 6.9 bar)
Storage Temp High	38 °C (100 °F)
Storage Temp Low	-18 °C (0 °F)

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/valves/53175-01>

Spec. Version B