

# Custom HPLC Column Packing Service - Preparative (21.2, 30 mm ID)

Specification Sheet (Part/REF # 79653)



## Product Specifications

|                         |                                                           |
|-------------------------|-----------------------------------------------------------|
| Stationary Phase        | Customer Supplied                                         |
| Format                  | HPLC Column                                               |
| Hardware Inner Diameter | Customer Specified<br>21.2 mm<br>30 mm                    |
| Hardware Length         | Customer Specified<br>50 mm<br>100 mm<br>150 mm<br>250 mm |
| Hardware Material       | Stainless Steel                                           |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79653>

Spec. Version B