

Semipreparative HPLC Column, [Custom Order] (7.0, 7.8, 10 mm ID)

Specification Sheet (Part/REF # 79642)



Inherent reversed-phase characteristics make the PRP-1 a rugged, multi-purpose column.

Product Specifications

Separation Mode	Reversed-Phase Anion Exchange Cation Exchange Ion Exclusion
Stationary Phase	(Choose a phase)
Format	HPLC Column
Hardware Inner Diameter	7.0 mm 7.8 mm 10 mm
Hardware Length	(Choose a length)
Particle Size	(Choose a particle size)
Hardware Material	Stainless Steel

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/laboratory-products/hplc-columns/79642>

Spec. Version B