

10 µL Conductive Filter Tips for ZEUS

Specification Sheet (Part/REF # 235901-ZEUS)



Hamilton offers a wide range of pipette tips to address a variety of common automation challenges.

Product Specifications

Sterilized	Non-Sterile
Tip Type	Standard
Tip Volume	10 µL
Maximum Tip Capacity	10 µL
Filter	Yes
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Light Red
Tip Material	Polypropylene - PP
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Rack Material	Polypropylene - PP
Tip Recycling Code	5
Rack Recycling Code	5
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
SLAS ANSI	N/A
Tip Length	29.9 mm
Orifice Size	0.3 mm

Case Size

Case of 5,760 tips (Blister 5 x 96 tips per rack)

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/oem/oem/235901-zeus>