

1.25 mL, Model 1002.5 TLL SYR, Instrument Syringe

Specification Sheet (Part/REF # 5495-40)



Product Specifications

Volume	1.25 mL
Syringe Type	Gastight Half-Height
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE, PFA
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.287 inches (7.29 mm)
Barrel Outer Diameter	0.375 inches (9.53 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/oem/syringes/5495-40>

Spec. Version B