

250 µL, Model 1750.5 TLLX SYR, Instrument Syringe

Specification Sheet (Part/REF # 5495-30)



Product Specifications

Volume	250 µL
Syringe Type	Gastight Half-Height
Syringe Series	1700
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE, PFA
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.128 inches (3.26 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/oem/syringes/5495-30>

Spec. Version B