

50 µL, Model 1710.5 TLLX SYR, Instrument Syringe

Specification Sheet (Part/REF # 5495-20)



Product Specifications

Volume	50 µL
Syringe Type	Gastight Half-Height
Syringe Series	1700
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, PTFE, PFA
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/oem/syringes/5495-20>

Spec. Version B