

# 1 mL, Model 1002.1 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 8250-37)



## Product Specifications

|                       |                                 |
|-----------------------|---------------------------------|
| Volume                | 1 mL                            |
| Syringe Type          | Longlife Gastight Half-Height   |
| Syringe Series        | 1000                            |
| Termination           | PTFE Luer Lock (TLL)            |
| Plunger Style         | X-Style                         |
| Fluid Path            | Borosilicate glass, UHMWPE, PFA |
| Autoclavable          | No                              |
| Gas Sterilizable      | Yes                             |
| Pressure              | 200 psig (13.8 bar)             |
| Maximum Temperature   | 80°C / 176°F                    |
| Minimum Temperature   | 10°C / 50°F                     |
| Barrel Inner Diameter | 0.287 inches (7.29 mm)          |
| Barrel Outer Diameter | 0.375 inches (9.53 mm)          |
| Sales Quantity        | 1 syringe per pack              |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/oem/syringes/8250-37>

Spec. Version B